









Conditioning instructions for Teknokroma capillary columns

Teknokroma capillary columns are preconditioned and individually tested. Each columns is provided with its own quality control test report.

- 1. Install the column and set the inlet carrier gas pressure based on the column dimensions (see table below). Use high purity carrier gas only (oxygen and moisture-free).
- 2. Check for leaks.

NOTE: At elevated temperatures traces of oxygen and humidity can damage irreversibly the column.

- 3. Purge the column with carrier gas at 40°C for at least 15 min. Raise the oven temperature up to the maximum isothermal temperature of the column (see column quality control report) or up to 10-30°C above the maximum operating temperature of the analysis -whichever is lower- at a rate of 2-5°C/min.
- 4. Hold this temperature for 1h-2h or until the baseline stabilizes.

NOTE: Working with highly sensitive detectors (ECD, MSD, etc.) might require longer conditioning periods.

NOTE: Columns with thick phase films may need an extended conditioning.





Typical Column Head Pressures*					
Length	I.D.				
meters	0.10 mm	0.20 mm	0.25 mm	0.32 mm	0.53 mm
10	25 psi	6 psi	4 psi	2 psi	0.9 psi
15	39 psi	9 psi	6 psi	4 psi	1.4 psi
25	67 psi	15 psi	9 psi	6 psi	2.3 psi
30	82 psi	18 psi	11 psi	7 psi	2.8 psi
50		32 psi	20 psi	12 psi	5 psi
60		39 psi	24 psi	14 psi	6 psi
105			40 psi	28 psi	9 psi

^{*}For optimal carrier gas linear velocity: 20 cm/s for He, and 40 cm/s for H2.

RELATED PRODUCTS

YOUR ANALYSIS MATTERS









SG-SGT-F0301 Triple Filter (O2/Moisture/Hydrocarbons)
SG-SGT-B0010-B8 1-Position Baseplate - 1/8" Brass

SG-SGT-CO1005 Click-On Inline Triple Trap (Oxygen/Moisture/Hydrocarbons) Cartidge-Stainless Steel Click-On Connectors 1/8" Brass