

AROCLOR 1254 IN ISOCTANE (10 PPM)

Column: **TRB-5MS**, P/N TR-520262

Size: 60m x 0.25mm x 0.25 μ m

Injection: splitless 60s, 270°C

Sample: 1 μ L Aroclors 1242, 1254 and 1260 standards in isooctane (10ppm), (1 ppm PCB30 and PCB209 internal standards)

Carrier Gas: He, 1mL/min


Program temperature: 70°C (1min) @ 30°C/min to 130°C @ 2.5°C/min to 300°C (15min)

Detector: MS KONIK-TECH, Modo EI+ (70 eV), SIM m/z 186, 222, 292, 326, 360, 394, 430, 464, 498 (50 ms), Source 140°C, Interface 300°C, Photomultiplier 1000V.

Chromatogram provided by José Antonio Muñoz from KONIK-TECH, S.A

Peak Name

- 1 PCB30 (IS)
- 2 PCB52
- 3 PCB44
- 4 PCB70
- 5 PCB66 / PCB95
- 6 PCB101
- 7 PCB99
- 8 PCB97
- 9 PCB87
- 10 PCB110
- 11 PCB149
- 12 PCB118
- 13 PCB153
- 14 PCB105
- 15 PCB138
- 16 PCB209 (IS)



TKG 1224

www.teknokroma.com

Teknokroma S. Coop. C. Ltda. - Camí de Can Calders, 14 - 08173 Sant Cugat del Vallès (Barcelona) Spain
Tel: +34) 936 748 800 Fax: +34) 936 752 405 e-mail: comercial@teknokroma.es