


SEPARATION OF HYDROCARBONS (FAST CHROMATOGRAPHY)

Column: **TRB-1**, P/N TR-110441
Dimensions: 10m x 0.10mm x 0.40µm
Injection: 0.5µL standard Hydrocarbons (0.95%/comp. in 2,2,4-Trimethylpentane), split 1:200, 200°C
Carrier gas: He, constant pressure 40 psi (275.6KPa).
Oven temperature: 190°C (Isothermal)
Detector: FID, 200°C

Chromatogram provided by J.I. Gómez Cívicos, M^ªA. Uguina Zamorano and J.L. Sotelo Sancho from Universidad Complutense de Madrid



TKG 1057