


PESTICIDES ANALYSIS

Column: Meta X5 P/N TR-820232
 Dimensions: 30m x 0.25mm x 0.25µm
 Injection: 1.0 µl standard, 10 ppm in Isooctane, splitless, 250°C
 Carrier gas: He, constant pressure, 9 psi (62 kPa)
 Oven program: 100°C (3,1min.) to 170°C @ 50°C/min. to 300°C(5,6min.)
 @ 5°C/min.
 Detector: MSD @ 280°C, scan 50-500 amus

Chromatogram supplied by J. Díaz from Chromatography Department, IQS.



TKG 1078

www.teknokroma.com