

PAHs Analysis (Benzo(j) Fluoranthene Resolved)

Column: TRB-50, P/N: TR-500232

Dimensions: 30m x 0.25mm x 0.25um

Injection: 300°C, pulsed splitless, pulse 20psi @ 0.3min, 40ml/min @ 0.2min


Liner: Gooseneck splitless

Carrier Gas: He, constant flow 1.6ml/min

Oven temperature: 90°C (1min) to 215°C (0.5min) @ 25°C/min to 235°C @ 4°C/min to 280°C @ 15°C/min to 320°C (20 min) @ 4°C/min

Detector: FID, 325°C

Sample: 1ul PAH mix, 20ug/ml each component



TKG 1282