

ALLERGENS IN COSMETICS

Column: **TRB-WAX**, P/N TR-140232

Size: 30m x 0.25mm x 0.25µm

Injection: 1 µl standard (25 µg/ml) in dichloromethane, splitless (60s), 250°C

Carrier Gas: He, 1mL/min

Program temperature: 32°C (5min) @ 4°C/min to 250°C (5min)


Detector: MS KONIK-TECH, Mode EI+ (70 eV), rango mas 35-300, Scan time 35ms,

Source Temperature 120°C, Interface 250°C, Photomultiplier 750V.

Chromatogram provided by Nieves Sarrion from KONIK-TECH, S.A (Barcelona).

Peak Name

- 1 α-Pinene
- 2 Camphene
- 3 β-Myrcene
- 4 α-Terpinene
- 5 Limonene
- 6 γ-Terpinene
- 7 Terpinolene
- 8 β-citronellal
- 9 Linalool
- 10 Linalyl acetate
- 11 Neral
- 12 s-Carvone
- 13 Geranial
- 14 β-Citronellol
- 15 α-
- 16 Geraniol
- 17 Benzyl alcohol
- 18 Hydroxycitronellal
- 19 Eugenol
- 20 Thymol
- 21 Cinnamyl alcohol
- 22 Farnesol (1)
- 23 Farnesol (2)
- 24 Isoeugenol
- 25 Hexyl cinnamal
- 26 Farnesol (3)
- 27 Coumarine
- 28 Benzyl benzoate
- IS Ethyl heptanoate



TKG 1199

www.teknokroma.com

Teknokroma S. Coop. C. Ltda. - Camí de Can Calders, 14 - 08173 Sant Cugat del Vallès (Barcelona) Spain
Tel: +34) 936 748 800 Fax: +34) 936 752 405 e-mail: comercial@teknokroma.es