

**TRB-BIODIESEL COLUMN**

Analysis of Glycerin and Glycerides (EN14105/ASTM D6584)  
Low Bleed at 370°C

Column: **TRB-BIODIESEL**, with retention gap, 2m x 0.53mm ID  
Injection: Biodiesel standard, cool on column  
Carrier gas: H<sub>2</sub>, constant flow 3 mL/min  
Program temperature: 50°C(1min) @ 15°C/min to 180°C @ 7°C/min to 230°C @ 30°C/min to 380°C(5min)  
Detector: FID, 380°C (N<sub>2</sub>, make up)

Chromatogram obtained from the analysis of Rapeseed oil. Provided by Mercedes H. Muñoz from ACCIONA

