

## COLUMNAS EQUIVALENTES USP

<b>CODIGO</b>	<b>DESCRIPCIÓN GENERAL</b>	<b>COLUMNA CAPILAR</b>
G1	Dimethylpolysiloxane oil	TRB-1, TRB-1ms
G2	Dimethylpolysiloxane gum	TRB-1, TRB-1ms
G3	50%phenyl-50%methylpolysiloxane	TRB-50
G5	3-cyanopropylsiloxane	TR-CN100
G8	90%-3-cyanopropyl-10%phenylmethylsiloxane	TR-CN100
G9	Methylvinylpolysiloxane	TRB-1, TRB-1ms
G14	Polyethylene glycol (MW 951-1050)	TR-WAX
G15	Polyethylene glycol (MW 3000-3070)	TR-WAX
G16	Polyethylene glycol (MW 15000)	TR-WAX
G19	25%phenyl-25%cyanopropylmethylsiloxane	TRB-225
G20	Polyethylene glycol (MW 380-420)	Meta.WAX 400
G25	Polyethylene glycol TPA	TRB-FFAP
G27	5%phenyl-95%methylpolysiloxane	TRB-5, TRB-5ms, Meta.X5
G28	25%phenyl-75%-dimethylpolysiloxane	TRB-20
G32	20%phenylmethyl-80%-dimethylpolysiloxane	TRB-20
G35	Polyethylene glycol with Nitroterephthalic acid	TRB-FFAP
G36	1%vinyl-5%phenylmethylpolysiloxane	TRB-5, TRB-5ms, Meta.X5
G39	Polyethylene glycol (MW=1500)	TR-WAX
G42	35%diphenyl-65%dimethylpolysiloxane	TRB-35
G43	6%cyanopropylphenyl-94%dimethylpolysiloxane	TRB-624, TRB-1301, TR-G43
G46	14%cyanopropylphenyl-86%dimethylpolysiloxane	TRB-1701