



## Removal systems for NS29/32 ground neck bottles

### SafetyCap ground neck

The locknut allows the stopper to be easily loosened and removed, even after long usage. For bottles with NS 29/32 mm standard ground joints.

- 1 to 2 connectors
- optional: SafeLock stopper → Pages 147 + 26

Part No.	Description	Connectors
107507	SafetyCap for NS29/32 ground neck bottles	without hole
107607	SafetyCap for NS29/32 ground neck bottles	2 x connector Ø 3.2 mm OD



GL45 adapter  
Part No. 107508



NS29/32 adapter  
Part No. 107509

## Adapter GL45 – NS29/32

Doesn't fit? Can't be. Continue to use original cap and simply adjust it to the bottle. With these adapters, it's fast and problem-free.

- optional: SafeLock stopper → Page 147

Part No.	Description
107508	GL45 adapter (f, bottle) on NS29/32 (f, cap)
107509	NS29/32 adapter (m, bottle) on GL45 (m, cap)

**About 12 million working hours are lost every year due to work-related illnesses**

