

A totally porous spherical packing, with a particle size of 4µm, giving a compromise alternative between the packings of 3 and 5µm.

## Novafix™ Cartridge SUPERSPHER

Function		Dimensions		
Particle size	(µm)	7.5 x 0.4 cm	15 x 0.4 cm	25 x 0.4 cm
Si 60	4	TR-015262	TR-015263	TR-015264
60 RP-8	4	TR-015265	TR-015266	TR-015267
100 RP-18	4	TR-015268	TR-015269	TR-015270
60 RP-8 ec	4	TR-015271	TR-015272	TR-015273
100 RP-18 ec	4	TR-015274	TR-015275	TR-015276

## Analytical columns SUPERSPHER

Function		Dimensions									
Particle size (µm)		10 x 0.46 cm	10 x 0.4 cm	12.5 x 0.46 cm	12.5 x 0.4 cm	15 x 0.46 cm	15 x 0.4 cm	20 x 0.46 cm	20 x 0.4 cm	25 x 0.46 cm	25 x 0.4 cm
Si 60	4	TR-011701	TR-411701	TR-011703	TR-411703	TR-011705	TR-411705	TR-011707	TR-411707	TR-011709	TR-411709
60 RP-8	4	TR-011711	TR-411711	TR-011713	TR-411713	TR-011715	TR-411715	TR-011717	TR-411717	TR-011719	TR-411719
100 RP-18	4	TR-011721	TR-411721	TR-011723	TR-411723	TR-011725	TR-411725	TR-011727	TR-411727	TR-011729	TR-411729
60 RP-8 ec	4	TR-011731	TR-411731	TR-011733	TR-411733	TR-011735	TR-411735	TR-011737	TR-411737	TR-011739	TR-411739
100 RP-18 ec	4	TR-011741	TR-411741	TR-011743	TR-411743	TR-011745	TR-411745	TR-011747	TR-411747	TR-011749	TR-411749