

Syringe Accessories Overview

Pipetting Aids/ Syringe Accessories

Reproducibility (Chaney) Adapter

The design allows reproducible dosing and also prevents the plunger from bending. The system can be retrofitted to 700, 800, 1700, 1800 and 7000 series syringes.

The Chaney adapter is available as an option for the above mentioned syringes and can be easily installed in the field.

Syringe Guide

The guide prevents the plunger from bending. Furthermore, it prevents the plunger from being pulled out. The guide is available for syringes series 700, 1700 and 7000. It is available separately and can be easily installed in the field.

Adapters

Septum Adapter

Store and sample syringes' contents with this adapter. It attaches to any TLL or LT syringes. The septum is replaceable.

RN Adapters (Removable Adapters)

There are three different removable adapters for 250 µl and larger syringes as well as for SampleLock Syringes available. It provides connection to threaded fittings, tubing assemblies, tees and crosses.

Syringe and Needle Cleaning Kit

The HAMILTON cleaning set helps to keep your syringe clean. Cleaned syringes last longer. The powerful cleaning solution and the detailed description of the cleaning procedure help to care your HAMILTON syringes.

GC Septa

These new long life GC Autosampler septa are designed to provide the longest life and lowest bleed possible when combined with HAMILTON GC Autosampler syringes.

MTB Sampling Valve

This lever-actuated on/off valve is used for storing and transporting samples in any GASTIGHT® TLL syringe.

Syringe Rack

Clean and safe storage of syringes from the series 700, 1700 and 7000 on your working place. Usable as free standing or wall mounted model.

