
















Tk Teknokroma Vials for Agilent

Autosampler Vials, Closures and Inserts for HPLC and GC for AGILENT

Features:

- Wide opening vial is easier to fill
- Fit most autosamplers that accommodate the 12x32 mm vial format
- Superior quality 33 expansion borosilicate clear glass (Type 1, Class A) or 51A amber glass
- Uniformly flat bottom for security with inserts
- I-D™ vials feature write-on patch with graduations at 0.5, 1.0 and 1.5 mL.
- Target Colorband™ crimp top vials feature color-coded graduated write-on patch
- Vials with an 11 mm crimp finish require aluminium seal closures
- Crimp/snap vials can be used with either 11 mm aluminium seals or Snap-It™ snap caps


Wide Opening Crimp Top Vials & Crimp Caps (12 x 32 mm):

Cat.No	P/N	Manufacturer	Description	Pk		
			CC-C4011-1	5181-3375	Vial 2 ml Clear, Wide Opening Crimp Top	100
			CC-C4011-1AP	5181-1210	11 mm Crimp Seal, with PTFE/Red Rubber, standar septum	100
			CC-C4011-1W	5182-0543	Vial 2 ml Clear, Wide Opening Crimp Top write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
			CC-C4011-2A	5181-1211	11 mm Crimp Seal, with red PTFE/White Silicone/Red PTFE, for aggressive solvents	100
			CC-C4011-2		Vial 2 ml Amber, Wide Opening Crimp Top	100
			CC-C4011-2W	5181-3376	Vial 2 ml Amber, Wide Opening Crimp Top write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
			TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100
			CC-C4011-LV1W		MacroVial 250 µl Fused Insert Vial, Clear Glass I-D, Write-on patch.	100
			TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100
			CC-03-FIV		300 µl Glass insert, fused into a 2 ml Crimp Top Vial, Write-on patch.	500
			TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100

Teknokroma Vials for Agilent **Tk**




Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-LV2W		MacroVial 250 ul. Fused Insert Vial, Amber Glass I-D, write-on patch	100
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100



CC-C4011-10		Solid Glass MicroVial 400 µl	12
TR-400012		11 mm Crimp Seal, with PTFE/Red Rubber standard	100

Crimp Seals With Different Septum Materials



CC-C4011-3A		11 mm Crimp Seal, with Polypropylene, for Corrosive Solvents	100
CC-C4011-4A		11 mm Crimp Seal, with PTFE/White Silicone, for High Purity Septum	100
CC-C4011-5A		11 mm Crimp Seal, with Black Viton Elastomer, for ECD, High Chemical Resistance	100
CC-C4011-1AP		11 mm Crimp Seal, with PTFE/Red Rubber, Standard Septum	100

Wide Opening Snap It & Snap Caps (12 x 32 mm):

Snap-it Vials can be used with Snap Caps or aluminium seals closures



Cat.No	P/N Manufacturer	Description	Pk
CC-C4011-5	5182-0544	Snap-It™ Vial 2 ml Clear, Wide Opening	100
CC-C4011-6		Snap-It™ Vial 2 ml Amber, Wide Opening	100



CC-C4011-5W	5182-0546	Snap-It™ Vial 2 ml Clear, Wide Opening. Write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
CC-C4011-51B	5182-3458	Blue Snap Cap with PTFE/Red Rubber	100



CC-C4011-6W	5182-0545	Snap-It™ Vial 2 ml Amber, Wide Opening. Write-on patch with graduations 0.5, 1.0 and 1.5 ml.	100
CC-C4011-54G	5182-3457	Green Snap Cap with red PTFE/white Silicone, High purity septum	100