




Zero Dead Volume Union **Tk**

High Pressure Union

PEEK ZDV Unions

Our PEEK zero-dead-volume (ZDV) Unions come complete with two F-300 Fingertight Fittings for 1/16" OD, tubing and are pressure rated to 5,000 psi (344 bar).




PEEK ZDV Unions

Cat.No.	Description
UP-P-742	PEEK ZDV, Union 0.010", thru-hole with F-300 Fittings
UP-P-704	PEEK ZDV, Union 0.020", thru-hole with F-300 Fittings
UP-P-760	PEEK ZDV, Union 0.050", thru-hole with F-300 Fittings

Stainless Steel ZDV Unions

Our high pressure, zero-dead-volume (ZDV) unions are precision machined from 316 stainless steel, carefully passivated, then thoroughly rinsed. Each comes complete with stainless steel nuts and ferrules.



Stainless Steel ZDV Unions

Cat.No.	Description
UP-U-402	ZDV, Union 0.020", thru-hole for 1/16" OD tubing
UP-U-435	ZDV, Union 0.010", thru-hole for 1/16" OD tubing
UP-U-438	ZDV, Union 0.067", thru-hole for 1/16" OD tubing
UP-U-411	NanoTight, Union 0.007", thru-hole for 1/16" OD tubing
UP-UH-402	VHP + Union 0.010", thru-hole for 1/16" OD tubing
UP-U-322	VICI Compatible ZDV Union 0.020" thru-hole for 1/16" OD
UP-U-412	Waters compatible ZDV Union 0.020" thru-hole for 1/16" OD