



Structure of Polyethylene glycol

TRB-WAXOmega Equivalent Phase

Supelco: Omegawax
Restek: Famewax

TRB-WAXOmega

Internal Length	Film	Temp	Part.	
Diam. (mm)	(m)	Thickness (µm)	limits (°C)	N°. (P/N)
0,25	30	0,25	40 to 260/270	TR-840232
0,32	30	0,25	40 to 260/270	TR-840233
0,53	30	0,50	40 to 260/270	TR-840535

TRB-WAXOmega

100% Polyethylene glycol, bonded and crosslinked phase.

- 100% Polyethylene glycol (PEG)
- High polarity column
- Specially designed for the analysis of Omega-3 and Omega-6 fatty acids methyl esters

