


GROB test

Column: TBR-5, 30m x 0.25 mm ID x 0.25 µm.
Carrier gas: He, 12 psi.
Oven: 110°C (Isothermal).
Injection: 1µl, split. (1:100), 260°C
Detector: FID, 280°C



TKG 1100

www.teknokroma.com

Teknokroma S. Coop. C. Ltda. - Camí de Can Calders, 14 - 08173 Sant Cugat del Vallès (Barcelona) Spain
Tel: +34) 936 748 800 Fax: +34) 936 752 405 e-mail: comercial@teknokroma.es